

Report Resolution of the PLANET



Introduction

The **planet** is the high resolution portable X-ray diffractometer from xplorex. With 'high resolution' we mean that the attainable resolution with the **planet** is comparable to a standard laboratory system and thus by far the best resolution hitherto obtained with a portable system.

In this note we demonstrate the attainable resolution with the **planet** with a normal measurement on normal commercial LaB₆.

Experimental

About 1g of the as-received powder was pressed into the sample holder (so-called front loading). The prepared specimen was mounted in the diffractometer and measured. Because the **planet** uses a reflection geometry, sample preparation is not critical even for heavy absorbers like LaB₆. The measurement parameters are summarized in Table 1.

| Parameter | Value |
|-----------------------------------|--|
| Range | |
| Start Angle ($^{\circ}2\theta$) | 19.3 |
| End Angle ($^{\circ}2\theta$) | 101.3 |
| Step size | Variable; The planet has predefined step sizes |
| Angle of incidence ($^{\circ}$) | 14.47 |
| Integration time | 120s/datapoint (the measurement comprises 14 times 640 datapoints) |
| Total Measurement time | 32 minutes |
| Diffractometer settings | |
| High Tension | 30 kV |
| Emission Current | 0.65 mA |
| Tube anode | Cu |
| Focus dimensions | 40µm diameter |
| Take off angle ($^{\circ}$) | 8 |
| Beam divergence ($^{\circ}$) | 1.25 |
| Specimen dimensions | 7mm diameter; 2mm thickness |
| Spinning frequency | 0.5 Hz |
| Optical path | Seemann – Bohlin based |
| Focusing circle radius | 160mm |
| Detector | Dectris' Mythen 1D solid state linear detector |
| Identification software | Match! From Crystal Impact |
| Reference database | Crystallographic Open Database |

Table 1 Measurement parameters for LaB₆ sample

The diffractometer settings are constant for the **planet**. The operator can choose the measurement range and the integration time for optimal results and ease-of-use.

Results

Figure 1 shows the complete diffraction pattern we recorded from the as received LaB₆ specimen.




Figure 1 Diffraction pattern from commercial LaB₆

To assess the attainable resolution of the **planet** we performed careful profile fitting on each observed reflection from LaB₆ and we reported the FWHM. The results of this exercise are shown in Table 2.

Table 2 Full Width at Half Maximum for LaB₆ reflections

| Scattering angle ($^{\circ}2\theta$) | FWHM ($^{\circ}2\theta$) |
|--|----------------------------|
| 21.357 | 0.198 |
| 30.384 | 0.109 |
| 37.441 | 0.092 |
| 43.506 | 0.084 |
| 48.956 | 0.070 |
| 53.987 | 0.079 |
| 63.217 | 0.083 |
| 67.546 | 0.084 |
| 71.744 | 0.091 |
| 75.842 | 0.089 |
| 83.843 | 0.102 |
| 87.789 | 0.087 |
| 99.639 | 0.083 |

As can be seen the FWHM is below 0.11° for all reflections except the first one (0.20°), where the beam path is significantly shorter. The length of the beam and the focus dimensions determine the attainable

resolution for any reflection. This confirms that the attainable resolution for all reflections is best in class, for any portable XRD in the market.

Comparison with a laboratory-based system

Because we used commercial LaB₆, we had to validate that this material is suitable for resolution measurements. To this end we measured both our commercial sample and a standard material (NIST 660a) on a laboratory-based system. These measurements are given by the Red (commercial) and Black (standard) curve in Figure 2.




Figure 2 FWHM measurements on laboratory system and planet

In the same graph we plotted the FWHM values we measured on the **planet** (Blue curve). We determined the FWHM by profile fitting in "Match!" from Crystal Impact. Please note that, although we measured the commercial sample with the **planet**, the FWHM values of all but the first reflections are even below the ones measured on the standard sample with a normal laboratory-based diffractometer! Moreover, we see the best resolution between 40 and 65°2θ, where it matters most, because generally the peak density is highest in that range.

Conclusion

The attainable resolution of the planet compares well to the attainable resolution of a standard laboratory-based instrument. The attainable resolution of the planet is best in class among currently available portable X-ray diffractometers.